



The physiological and biochemical mechanisms of action of V.A.C.[®] Therapy[™]

Presented by
Dr. Tino Hauser
Director Clinical Marketing

Overview

- Summarize experimental evidence for the mechanism of action of V.A.C.[®] Therapy[™]

- Published literature
 - Theoretical models
 - In vitro studies
 - Animal experiments
 - Clinical studies

- Key conclusions regarding the role of V.A.C.[®] Therapy[™] in the wound healing process

Chronic versus acute wound healing

Healing wounds



Mitogenic activity

Inflammatory Cytokines

Protease activity


Mitotically competent cells

Chronic wounds




Teot L, Banwell PE, Ziegler UE, Surgery in Wounds. Springer 2004, 30


Possible mechanisms of action of V.A.C.[®] Therapy[™]




Possible mechanisms of action of V.A.C.[®] Therapy[™]




V.A.C.® Therapy™ – Foam/Tissue Interactions, Forces on Tissue Surface




Therapy initiated:



Equilibrium:



Foam/Tissue interactions and Wound Tissue Stretching




- ➡ Negative Pressure (Vacuum)
- ➡ Foam **Normal** Compressive Stress
- ➡ Foam **Lateral** Compressive Stress

V.A.C.® Therapy™ induces microdeformations

- Microdeformations in V.A.C.® treated wounds
- V.A.C.® treated wounds show a rich vascular network compared to control sites without foam




Non-V.A.C.® System (day 7)




V.A.C.® System (day 7)

Saxena V, et al. *Plastic and Reconstructive Surgery*, 2004; 114(5):1086-96.

V.A.C.® Therapy™ induces microdeformations




V.A.C.® System
↓
Microdeformation



Common suction device
↓
No microdeformation


Orgill DP, Bayer LR, Neuwalder J, Felter RC. Business Briefing: Global Surgery – Future Directions 2005, 22-25

V.A.C.® Therapy™ induces wound surface strain in a repeating pattern




Modified Graphic from Willy C, *Die Vakuumtherapie*, May 2005, p86

V.A.C.[®] Therapy[™] may directly stimulate cell proliferation by microdeformations



Saxena V, et al. Plastic and Reconstructive Surgery, 2004; 114(5):1086-96. Huang S, Chen CS, Ingber DE. Mol Biol Cell. 1998 Nov;9(11):3179-93. Danciu TE, Gagari E, Adam RM, Damoulis PD, Freeman MR. J Dent Res. 2004 Aug;83(8):596-601. Katsumi A, Orr AW, Tzima E, Schwartz MA. J Biol Chem. 2004 Mar 26;279(13):12001-4. Epub 2004 Feb 11.

V.A.C.® Therapy™ may directly stimulate cell proliferation



Microvessel density of wounds treated with V.A.C.® Therapy™ and GranuFoam® contact was significantly higher compared to areas not covered with foam during 1st and 2nd week of treatment. Wounds treated with V.A.C.® Therapy™ and GranuFoam® had greater microvessel density compared with the same wound prior to treatment (p=0.02).

Greene AK, Ann Plast Surg. 2006 Apr;56(4):418-422.

V.A.C.® Therapy™ may directly stimulate cell proliferation

V.A.C.® Therapy™ activates key players of cell cycle regulation

Method:

In vitro study with cultured fibroblasts


Cells were stretched (5% for 24h, 48h, and 72h) with 5 min stretch/2 min relax cycle

Result:

a) Phosphorylated **c-Jun** ↑ (at 48 and 72h of stretching)


b) Phosphorylated **p38** ↑ (at 24, 48 and 72h of stretching)

c-Jun belongs to JNK signaling pathway; p38 belongs to MAPK signaling pathway




Kremers L, Kearns M, Hammon D, Scott AC, Daniel L, Morykwas MJ. Wound Repair and Regeneration, 11(5),A3.


Possible mechanisms of action of V.A.C.® Therapy™



Role of local perfusion in wound healing




V.A.C.® Therapy™ increases local blood perfusion (humans/healthy skin)



- A. In intact human forearm skin mean cutaneous blood flow under the foam increases significantly when V.A.C.® Therapy™ negative pressure increases to up to -300 mmHg: $p < 0.001$
- B. At -300 mmHg: over 5-fold increase in perfusion with GranuFoam® and nearly 3-fold increase with VersFoam®
- C. No statistically significant increase was observed when pressure was raised above 300 mmHg

Timmers et al. Ann Plast Surg. 2005 Dec;55(6):665-71.


V.A.C.® Therapy™ increases local blood perfusion (animals/ wounds)



- Peak flow of four times baseline values at -125 mmHg
- Blood flow was depressed below baseline at negative pressures of -400 mmHg and above
- Intermittent application resulted in repeated increases in blood flow

Morykwas et al. Ann Plast Surg. 1997 Jun;38(6):553-62.

V.A.C.[®] Therapy[™] increases local blood perfusion (in human burn wounds)



Seven patients with bilateral partial thickness hand burns were treated with V.A.C.[®] Therapy[™] on the more intense injured hand and with silver sulphadiazine cream on the contra lateral one

On day 2 and 3 perfusion was significantly better in the V.A.C.[®] Therapy[™] treated hand (resp. 20.1% (p = 0.001) and 15.2 % (p = 0.006). Oedema was reduced in the V.A.C.[®] Therapy[™] treated hand

Kamolz et al. Burns 2004 (30):253-58.


V.A.C.® Therapy™ reduces localized edema

- Edema reduction in hand burns
- Massive reduction of edema in the V.A.C.® Therapy™ group (picture left hand)
- Very large amount of fluid (up to 500 ml) was removed with V.A.C.® Therapy™
- Potential Mechanisms: *Mechanical removal of fluid and effects on the inflammatory response that causes edema*
- Clinical Effect: *Reduced tissue pressure and improved capillary blood flow*

Kamolz LP et al., Burns. 2004 May;30(3):253-8.




Possible mechanisms of action of V.A.C.® Therapy™




Role of granulation tissue formation in wound healing

- The formation of granulation tissue in an open wound allows the re-epithelialization phase to take place, as epithelial cells migrate across the new tissue to form a barrier between the wound and the environment.
- Granulation tissue is thus needed to fill the void that has been left by a large, open wound that crosses the basal membrane.
- Absence of granulation tissue is a feature of a non-healing wound.



Reviewed by Romo and McLaughlin, <http://www.emedicine.com/ent/topic13.htm> 2003. Lazarus et al., Arch Dermatol 1994;130:489-93.


V.A.C.[®] Therapy[™] promotes granulation tissue formation



- A. V.A.C.[®] Therapy[™] treated wounds filled with granulation tissue at a significantly greater rate than control wounds (with both continuous and intermittent therapy), $p \leq 0.01$
- B. The mean increase with continuous negative pressure (125 mmHg) was 63.3%, the mean increase for intermittently treated wounds was 103.4%

Morykwas et al. Ann Plast Surg. 1997 Jun;38(6):553-62.

V.A.C.® Therapy™ promotes granulation tissue formation




Highest wound volume reduction and rate of granulation tissue formation in wounds treated with standard recommended negative pressure of -125 mmHg. Wounds treated with V.A.C.® Therapy™ (-125 mmHg) filled with a significant higher rate than wounds treated with low (-25 mmHg) or high (-500 mmHg) negative pressure ($p < 0.0001$)

Unregulated air leakage leads to wound volume increase and is a serious risk

*Wounds were treated until one of the wounds has granulated to a level flush with the surrounding tissue. Wounds treated with 125mm Hg had filled 1st completely with granulation tissue (by day 8).

Morykwas et al. Ann Plast Surg. 2001 Nov;47(5):547-51.

Possible mechanisms of action of V.A.C.® Therapy™



V.A.C.® Therapy™ removes wound exudate containing molecules that can inhibit wound healing

The Role of Matrix Metalloproteinases in Wound Healing

- Degradation of all major components of extracellular matrix
- Important for normal wound healing (e.g. removal of damaged tissue cell migration)
- MMP family
 - Collagenases (MMP-1, MMP-8)
 - Gelatinases (MMP-2, MMP-9)
 - Strong collagenase and gelatinase (MMP-13)

Higher levels in chronic wound fluid (versus acute wounds)



Degradation of key molecules for wound healing

The Role of proinflammatory Cytokines in Wound Healing

- Non-antibody proteins released by one cell population (contact with specific antigen)
- Intercellular mediators for e.g. immune response
- TNF- α is an important and well investigated proinflammatory cytokine

Higher levels of TNF- α in chronic wound fluid (versus acute wounds)




- Higher MMP activity
- Inhibition of Cell Cycle

Murphy G et al. Biochem Soc Trans. 1990 Oct;18(5):812-5. Agren MS et al. J Invest Dermatol. 1992 Dec;99(6):709-14. Teot L, Banwell PE, Ziegler UE, Surgery in Wounds. Springer 2004, 30

V.A.C.® Therapy™ reduces activity of selected MMPs and inflammatory cytokines

Effect of V.A.C.® Therapy™ with and without GranuFoam® contact on MMP-2/9 activity¹



Role of MMP-2 and MMP-9²

- 5-10 fold increased in chronic versus acute wound fluid
- Matrix degradation and hypothetical effect on endothelial proliferation (lowering MMP-mediated angiostatin and endostatin production)

1. Greene AK et al. Ann Plast Surg. 2006 Apr;56(4):418-422 2. Wysocki AB et al. J Invest Dermatol. 1993 Jul;101(1):64-8.

V.A.C.® Therapy™ reduces activity of selected MMPs and inflammatory cytokines


Effect of V.A.C.® Therapy™ on the Expression of MMP-1, MMP-2 and MMP13 in Human Granulation wounds

- Clinical study with n = 5 patients (chronic wounds)
- MMP-1/2/13 expression decreased following V.A.C.® Therapy™ treatment (1 week after treatment start)
- MMP-13 showed most significant decrease ($p < 0.05$), mRNA levels were reduced to 1/10 at day 4 and 7



Shi B et al. Zhonghua Zheng Xing Wai Ke Za Zhi. 2003 Jul;19(4):279-81

V.A.C.® Therapy™ reduces activity of selected MMPs and inflammatory cytokines




V.A.C.® Therapy™ Reduces TNFα and local Inflammation¹

- **Goal:** Determine if V.A.C.® Therapy™ has an effect on inflammatory cytokines, proteases and protease inhibitors
- **Result:** Hospitalized pressure ulcer patients (n=8) have a significant decrease in tumor necrosis factor-alpha (TNFα) following application of V.A.C.® Therapy™

1. Stechmiller et al., to be published in May – June, 2006, *Wound Repair & Regeneration*

Possible mechanisms of action of V.A.C.® Therapy™



V.A.C.® Therapy™ provides a closed moist wound environment

- Generally accepted that moist wound environment supports healing
- Adhesive semi-occlusive drape retains moisture and allows gas exchange
- Prevents tissue dehydration
- Prevents cell death
- Adhesive semi-occlusive drape protects wound and reduces the chance of wound contamination

Summary

- Framework of basic mechanisms is supported by a theoretical model and experimental evidence
- V.A.C.® Therapy™ supports wound healing by several basic mechanisms
- Different mechanisms of action are interrelated and can act synergistically
- Major clinical effects (faster wound healing) can be explained by the evidence-based framework of basic mechanisms

References

1. Agren MS, Taplin CJ, Woessner JF Jr, Eaglstein WH, Mertz PM. Collagenase in wound healing: effect of wound age and type. *J Invest Dermatol*. 1992 Dec;99(6):709-14.
2. Chen C, Schultz GS, Bloch M, Edwards PD, Tubes S, Mast BA. Molecular and mechanistic validation of delayed healing rat wounds as a model for human chronic wounds. *Wound Repair Regen*. 1999 Nov-Dec;7(6):486-94.
3. Danciu TE, Gagari E, Adam RM, Damoullis PD, Freeman MR. Mechanical strain delivers anti-apoptotic and proliferative signals to gingival fibroblasts. *J Dent Res*. 2004 Aug;83(8):596-601.
4. Edwards, Ruth; Harding, Keith G Bacteria and wound healing. *Current Opinion in Infectious Diseases*. 17(2):91-96, April 2004.
5. Ford CN, Reinhard ER, Yeh D, Syrek D, De Las Morenas A, Bergman SB, Williams S, Hamori CA. Interim analysis of a prospective, randomized trial of vacuum-assisted closure versus the Healthpoint system in the management of pressure ulcers. *Annals of Plastic Surgery*. 2002 Jul;49(1):55-61; discussion 61.
6. Gladden, L. B. Lactate metabolism: a new paradigm for the third millennium. *J Physiol*, 558: 5-30, 2004.
7. Greene AK, Puder M, Roy R, Arsenault D, Kwei S, Moses MA, Orgill DP. Microdeformational Wound Therapy: Effects on Angiogenesis and Matrix Metalloproteinases in Chronic Wounds of 3 Debilitated Patients. *Ann Plast Surg*. 2006 Apr;56(4):418-422.
8. Hopf, H. W., Gibson, J. J., Angeles, A. P., Constant, J. S., Feng, J. J., Rollins, M. D., Zamirul, H. M., and Hunt, T. K. Hyperoxia and angiogenesis. *Wound Repair Regen*, 13: 558-564, 2005.
9. Huang S, Chen CS, Ingber DE. Control of cyclin D1, p27(Kip1), and cell cycle progression in human capillary endothelial cells by cell shape and cytoskeletal tension. *Mol Biol Cell*. 1998 Nov;9(11):3179-93.
10. Kamolz LP, Andel H, Haslik W, Winter W, Meissl G, Frey M. Use of subatmospheric pressure therapy to prevent burn wound progression in human: first experiences. *Burns*. 2004 May;30(3):253-8.
11. Katsumi A, Orr AW, Tzima E, Schwartz MA. Integrins in mechanotransduction. *J Biol Chem*. 2004 Mar 26;279(13):12001-4. Epub 2004 Feb 11.
12. Kilpadi et al., to be published in March – April 2006, *Wound Repair and Regeneration*.

References

13. Kremers L, Kearns M, Hammon D, Scott AC, Daniel L, Morykwas MJ. Involvement of mitogen activated protein kinases (MAP kinases) in increased wound healing during sub-atmospheric pressure (SAP) treatment. *Wound Repair and Regeneration*, 11(5),A3.
14. Lazarus, G. S., Cooper, D. M., Knighton, D. R., Margolis, D. J., Pecoraro, R. E., Rodeheaver, G., and Robson, M. C. Definitions and guidelines for assessment of wounds and evaluation of healing. *Arch Dermatol*, 130: 489-493, 1994.
15. Morykwas, M. J., Argenta, L. C., Shelton-Brown, E. I., and McGuirt, W. Vacuum-assisted closure: a new method for wound control and treatment: animal studies and basic foundation. *Ann Plast Surg*, 38: 553-562, 1997.
16. Morykwas, M. J., Faler, B. J., Pearce, D. J., and Argenta, L. C. Effects of varying levels of subatmospheric pressure on the rate of granulation tissue formation in experimental wounds in swine. *Ann Plast Surg*, 47: 547-551, 2001.
17. Murphy G, Hembry RM, Hughes CE, Fosang AJ, Hardingham TE. Role and regulation of metalloproteinases in connective tissue turnover. *Biochem Soc Trans*. 1990 Oct;18(5):812-5.
18. Orgill DP, Bayer LR, Neuwalder J, Felter RC. Microdeformational Wound Therapy – A New Era in Wound healing. *Business Briefing: Global Surgery – Future Directions* 2005, 22-25.
19. Payne, w. g., Wright, T. E., Ko, F., Wheeler, C., Wang, X., and Robson, M. C. Bacterial Degradation of Growth Factors. *the Journal of Applied Research in Clinical and Experimental Therapeutics*, 3: 35-40, 2006.
20. Robson, M. C. Wound infection. A failure of wound healing caused by an imbalance of bacteria. *Surg Clin North Am*, 77: 637-650, 1997.
21. Robson, M. C., Stenberg, B. D., and Heggens, J. P. Wound healing alterations caused by infection. *Clin Plast Surg*, 17: 485-492, 1990.
22. Romo, T. I. and McLaughlin, L. A. Wound healing, skin. <http://www.emedicine.com/ent/topic13.htm>, 2003, 2006.

References

23. Saxena V, Hwang CW, Huang S, Eichbaum Q, Ingber D, Orgill DP. Vacuum-assisted closure: microdeformations of wounds and cell proliferation. *Plast Reconstr Surg*. 2004 Oct;114(5):1086-96; discussion 1097-8
24. Shi B, Chen SZ, Zhang P, Li JQ. Effects of vacuum-assisted closure (VAC) on the expressions of MMP-1, 2, 13 in human granulation wound. *Zhonghua Zheng Xing Wai Ke Za Zhi*. 2003 Jul;19(4):279-81
25. Stechmiller et al., to be published in May – June, 2006, *Wound Repair & Regeneration*
26. Teot L, Banwell PE, Ziegler UE, *Surgery in Wounds*. Springer 2004, 30
27. Timmers MS, Le Cessie S, Banwell P, Jukema GN. The effects of varying degrees of pressure delivered by negative-pressure wound therapy on skin perfusion. *Ann Plast Surg*. 2005 Dec;55(6):665-71.
28. Wackenfors, A., Gustafsson, R., Sjogren, J., Algotsson, L., Ingemansson, R., and Malmso, M. Blood flow responses in the peristernal thoracic wall during vacuum-assisted closure therapy. *Ann Thorac Surg*, 79: 1724-1730, 2005.
29. Willy C, *Die Vakuumtherapie*, May 2005, ISBN 3-00-016219-4
30. Wysocki AB, Staiano-Coico L, Grinnell F. Wound fluid from chronic leg ulcers contains elevated levels of metalloproteinases MMP-2 and MMP-9. *J Invest Dermatol*. 1993 Jul;101(1):64-8.